

# Preparation of Competent Cells

(Reference: Inoue, Hiroaki, Nojima, Hiroshi, and Okayama, Hiroto. "High Efficiency Transformation of Escherichia coli with Plasmids", Gene, 96, 23-28, 1990)

1. Streak frozen stock of SURE cells on an LB + Tet agar plate.
2. Allow the cells to grow overnight at 37 °C
3. Pick 10 large (~ 2-3 mm in diameter) colonies and place all 10 colonies into 5 ml SOB.
4. Add the SOB medium containing the colonies to 250 ml of SOB medium in a 2-l flask.  
(Note: 2.5 ml of the magnesium stock solution should be added to the SOB medium in the 2-l flask prior to this step)
5. Incubate the flask at room temperature with vigorous shaking.
6. Grow cells to  $A_{600} = 0.6$ ; the specific growth rate,  $m$ , of SURE cells at R.T.  $\sim 0.44 \text{ h}^{-1}$ , and the double time at R.T. is about 1.6 h. It will take about 24 h to reach an O.D. of 0.6, but it will depend on the sizes of the colonies you picked . If you inoculate the cells say at 2 p.m., you should check the  $A_{600}$  in the morning when you come in, then you can calculate the exact time the cells will reach an O.D. of 0.6 using the following formula:  $\text{Time (h)} = 1/m \cdot \ln (0.6/A_{600})$
7. Place flask on ice for 10 min.
8. Transfer culture into a cooled centrifuge bottle and spin at 3500 rpm ( $\sim 2100 \times g$ ) in the JA10 rotor for 10 min at 4 °C.
9. Remove the supernatant and resuspend the pelleted cells in 80 ml of ice-cold buffer TB with a glass pipette; incubate on ice for 10 min.
10. Spin the cells down again in the same bottle at 3500 rpm for 10 min at 4°C.
11. Remove the supernatant and resuspend the pelleted cells in 20 ml of ice-cold buffer TB.  
Make sure the cells are well resuspended.
12. Add DMSO to cell suspension to 7% (i.e. 1.4 ml).
13. Incubate the cells on ice for another 10 min.
14. Aliquot the cell suspension into 250  $\mu\text{l}$  per eppendorff tubes and freeze the tubes quickly with liquid nitrogen.

SOB Medium:

(2% (w/v) Bacto Tryptone
0.5% (w/v) Yeast extract
10 mM NaCl

2.5 mM KCl
10 mM MgCl <sub>2</sub>
10 mM MgSO <sub>4</sub> )

To make 500 ml:

Bacto Tryptone	10 g
Yeast Extract	2.5 g
NaCl	0.29 g
KCl	0.09 g

- dissolve in 500 ml of H<sub>2</sub>O
- autoclave

Magnesium Stock Solution

1 M MgCl <sub>2</sub>	2 g
1 M MgSO <sub>4</sub>	2.5 g

- dissolve in 10 ml H<sub>2</sub>O
- filter sterilize

Before use, add 1/100 vol of Magnesium stock solution to medium (i.e. 2.5 ml to 250 ml medium).

Buffer TB:

10 mM Hepes  
15 mM CaCl<sub>2</sub>  
250 mM KCl  
55 mM MnCl<sub>2</sub>

To make 100 ml:

Hepes (M.W. 238.3)	0.24 g
CaCl <sub>2</sub> (M.W. 147.02)	0.22 g
KCl (M.W. 74.55)	1.86 g

- dissolve in 80 ml of H<sub>2</sub>O

- adjust pH to 6.7 with KOH
- then add 1.09 g of MnCl<sub>2</sub> (M.W. 197.9) and adjust the final volume to 100 ml with H<sub>2</sub>O
- filter sterilize the solution