

## MEMBRANE-FREE RIBOSOMES (mfrs)

The removal of endogenous vesicles also removes a large fraction of ribosomes, which must be added back to the reaction exogenously. The procedure for isolation of mfrs is as follows:

1. Inoculate 1.5L of TB + 0.2% glucose with 100ml saturated culture of MRE 600 cells.
2. Grow to OD600 of 1.0 at 30 C.
3. Pellet cells (5000 rpm, 15 minutes), and wash 1X with ddH<sub>2</sub>O.
4. Resuspend in buffer A at 0.5g/ml
5. Pass suspension through French pressure cell 3X at 4000psi
6. Centrifuge 30 min. at 40 000xg (19 000 rpm in 50.2).
7. Dilute supernatant 10- fold in buffer B.
8. Centrifuge 2 hours at 150 000xg (37 000 rpm in Ti 50.2).
9. Resuspend pellet in 18 mL of buffer B (Note: use Dounce homogenizer with type A pestle). Incubate overnight with rotation at 4 C.
10. Centrifuge 2 hours at 150 000xg.
11. Resuspend in 20 mL of buffer C.
12. Centrifuge 2 hours at 150 000xg. Note: Pellet is almost clean, and has approximate consistency of "snot".
13. Resuspend pellet in 500 $\mu$ l buffer C.

**Note:** Volume of mfrs required must be determined, but should be optimal in this case at about 0.1 $\mu$ l per 10 $\mu$ l reaction.